

SureSelect Target Enrichment System

Illumina Single-End Sequencing Platform Library Prep

Protocol

Version 1.2, April 2009

**SureSelect platform manufactured with Agilent
SurePrint Technology**

**Research Use Only. Not for use in Diagnostic
Procedures.**



Agilent Technologies

Notices

© Agilent Technologies, Inc. 2009

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

G3360-90010

Edition

Version 1.2, April 2009

Printed in USA

Agilent Technologies, Inc.
5301 Stevens Creek Rd
Santa Clara, CA 95051 USA

Acknowledgement

Oligonucleotide sequences © 2006 and 2008 Illumina, Inc. All rights reserved. Only for use with the Illumina Genome Analyzer and associated assays.

Technical Support

Technical product support may be obtained by contacting your local Agilent Support Services representative. Agilent's worldwide sales and support center telephone numbers can be obtained at the following Web site:

www.agilent.com/chem/contactus

or send an email to:

lsca-ibs-support@agilent.com

Notice to Purchaser

Research Use Only. Not for use in diagnostic procedures.

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

U.S. Government Restricted Rights. Software and technical data rights granted to the federal government include only those rights customarily provided to end user customers. Agilent provides this customary commercial license in Software and technical data pursuant to FAR 12.211 (Technical Data) and 12.212 (Computer Software) and, for the Department of Defense, DFARS 252.227-7015 (Technical Data - Commercial Items) and DFARS 227.7202-3 (Rights in Commercial Computer Software or Computer Software Documentation).

Safety Notices

CAUTION

A **CAUTION** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A **WARNING** notice denotes a hazard. It calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood and met.

In this Guide...

This guide describes Agilent's recommended operational procedures to capture genomic regions of interest using Agilent's SureSelect Target Enrichment System Kit and sample preparation kits for next-generation sequencing. This protocol is specifically developed and optimized to use Biotinylated RNA oligomer libraries, or Bait, to enrich targeted regions of the genome from repetitive sequences and sequences unrelated to the research focus.

This guide uses the Illumina single-end sequencing platform for library preparation.

If you have comments about this protocol, send an e-mail to feedback_genomics@agilent.com.

1 Before You Begin

This chapter contains information (such as procedural notes, safety information, required reagents and equipment) that you should read and understand before you start an experiment.

3 Hybridization

This chapter describes the steps to prepare and hybridize samples.

4 Post-Hybridization Amplification


This chapter describes the steps to amplify, purify, and assess quality of the sample library.

Content

1	Before You Begin	7
	Procedural Notes	8
	Safety Notes	8
	Required Reagents	9
	Optional Reagents	10
	Required Equipment	11
	Optional Equipment	12
2	Sample Preparation	13
	Step 1. Shear DNA	16
	Step 2. Purify the sample with the QIAquick PCR Purification Kit	17
	Step 3. Assess quality with Agilent 2100 Bioanalyzer	18
	Step 4. Repair the ends	19
	Step 5. Purify the repaired DNA with the QIAquick PCR Purification Kit	20
	Step 6. Add 'A' Bases to the 3' end of the DNA fragments	21
	Step 7. Purify the sample with a Qiagen MinElute PCR Purification Column	22
	Step 8. Ligate the adapter	23
	Step 9. Purify the sample with the QIAquick PCR Purification Kit	24
	Step 10. Select size using a gel	25
	Step 11. Purify the gel	26
	Step 12. Amplify prepped library	28
	Step 13. Purify the sample	29
	Step 14. Check gDNA quantity and quality	30
	Step 15. Assess quality with Agilent 2100 Bioanalyzer	31
3	Hybridization	33
	Step 1. Hybridize the library	36
	Step 2. Prepare magnetic beads	40
	Step 3. Select hybrid capture with SureSelect	41
	Step 4. Desalt capture solution	42

Contents

4	Post-Hybridization Amplification	43
	Step 1. Amplify the sample	44
	Step 2. Purify PCR	46
	Step 3. Assess quality with Agilent 2100 Bioanalyzer	47
5	Reference	49
	Alternative Capture Equipment Combinations	50



1 Before You Begin

Procedural Notes	8
Safety Notes	8
Required Reagents	9
Optional Reagents	10
Required Equipment	11

Make sure you read and understand the information in this chapter and have the necessary equipment and reagents listed before you start an experiment.

NOTE

Agilent cannot guarantee the Target Enrichment Kit and cannot provide technical support for the use of non-Agilent protocols to process samples for enrichment.



Procedural Notes

- To prevent contamination of reagents by nucleases, always wear powder-free laboratory gloves, and use dedicated solutions and pipettors with nuclease-free aerosol-resistant tips.
- Maintain a clean work area.
- Do not mix stock solutions and reactions containing gDNA on a vortex mixer. Instead, gently tap the tube with your finger to mix the sample.
- Avoid repeated freeze-thaw cycles of stock and diluted gDNA solutions.
- When preparing frozen reagent stock solutions for use:
 - 1 Thaw the aliquot as rapidly as possible without heating above room temperature.
 - 2 Mix briefly on a vortex mixer, then spin in a centrifuge for 5 to 10 seconds to drive the contents off of walls and lid.
 - 3 Store on ice or in a cold block until use.
- In general, follow Biosafety Level 1 (BL1) safety rules.

Safety Notes

CAUTION

- Wear appropriate personal protective equipment (PPE) when working in the laboratory.
-

Required Reagents

Table 1 Required Reagents for Library Prep and Post-Hybridization Amplification

Description	Vendor and part number
Agilent DNA 1000 Kit	Agilent p/n 5067-1504
Nuclease-free Water (not DEPC-treated)	Ambion, Cat # AM9930
MinElute PCR Purification Kit	Qiagen, Cat # 28004
Illumina Genomic DNA Sample Prep Kit	Illumina, Cat # FC-102-1001
Sybr-Gold	Invitrogen p/n S-11494
Trackit Cyan/Orange Loading Dye (50 mM Tris pH 8.0, 40 mM EDTA, 40% (w/v) sucrose)	Invitrogen p/n 10482-028
Trackit 50 BP DNA ladder	Invitrogen p/n 10488-043
NuSeive GTG Agarose	Lonza p/n 50080
QIAquick PCR Purification Kit	Qiagen, Cat # 28104
QIAquick or MinElute Gel Extraction Kit	Qiagen p/n 28704 or Qiagen p/n 28604
Herculase II Fusion DNA Polymerase	Stratagene p/n 600677
50x TAE buffer	
Distilled water	
100% Ethanol, molecular biology grade	
Isopropanol	

1 Before You Begin

Optional Reagents

Table 2 Required Reagents for Hybridization

Description	Vendor and part number
SureSelect Target Enrichment, 10 reactions	Agilent p/n G3360A
SureSelect Target Enrichment, 25 reactions	Agilent p/n G3360B
SureSelect Target Enrichment, 50 reactions	Agilent p/n G3360C
SureSelect Target Enrichment, 100 reactions	Agilent p/n G3360D
SureSelect Target Enrichment, 250 reactions	Agilent p/n G3360E
SureSelect Target Enrichment, 500 reactions	Agilent p/n G3360F
SureSelect Target Enrichment, 1000 reactions	Agilent p/n G3360G
SureSelect Target Enrichment, 2000 reactions	Agilent p/n G3360H
SureSelect Target Enrichment, 5000 reactions	Agilent p/n G3360J
Nuclease-free Water (not DEPC-treated)	Ambion, Cat # AM9930
Dynabeads [®] M-280 Streptavidin	Invitrogen, Cat # 112-05D

Optional Reagents

Table 3 Optional Reagents

Description	Vendor and part number
SureSelect Target Enrichment, Demo	Agilent p/n G4459A

Required Equipment

Table 4 Required Equipment for Library Prep and Post-Hybridization Amplification

Description	Vendor and part number
Agilent 2100 Bioanalyzer	Agilent p/n G2938C
Nuclease-free 1.5 mL microfuge tubes (sustainable at 95°C)	Ambion p/n AM12400 or equivalent
Dark Reader transilluminator	Clare Chemical Research, Inc. p/n DR45M
Covaris S-series Single Tube Sample Preparation System, Model S2	Covaris
Covaris microTUBE with AFA fiber and snap cap	Covaris p/n 520045
Nuclease-free 0.2 mL PCR tubes, thin-walled	Eppendorf p/n 951010006 or equivalent
Microcentrifuge	Eppendorf Microcentrifuge Model 5417C
UV-VIS spectrophotometer	NanoDrop p/n ND-1000 or equivalent
P10, P20, P200 and P1000 pipettes	Pipetman P10, P20, P200, P1000 or equivalent
Stratagene Gradient Cycler	Stratagene p/n G5100B
Disposable scalpels or razor blades	
Electrophoresis unit	
Electrophoresis power supply	
Gel trays and tank	
Ice bucket	
Powder-free gloves	
Sterile, nuclease-free aerosol barrier pipette tips	
Timer	
Vortex mixer	

1 Before You Begin

Optional Equipment


Table 5 Required Equipment for Hybridization

Description	Vendor and part number
Mx3000P/Mx3005P 96-well tube plates	Agilent p/n 410088
Mx3000P/Mx3005P optical strip caps	Agilent p/n 401425
BD Clay Adams Nutator Mixer	BD Diagnostics p/n 421105 or equivalent
Dynal DynaMag-2	Invitrogen p/n 123-21D or equivalent
P10, P20, P200 and P1000 pipettes	Pipetman P10, P20, P200, P1000 or equivalent
Pipet-Light Multichannel Pipette, 12 channels	Rainin p/n L12-20 or equivalent
Stratagene Gradient Cycler	Stratagene p/n G5100B
Sterile, nuclease-free aerosol barrier pipette tips	
Timer	
Vortex mixer	
Water bath	

Optional Equipment

Table 6 Optional for Equipment Hybridization

Description	Vendor and part number
Tube-stirp capping tool	Agilent p/n 410099



2 Sample Preparation

- Step 1. Shear DNA 16
- Step 2. Purify the sample with the QIAquick PCR Purification Kit 17
- Step 3. Assess quality with Agilent 2100 Bioanalyzer 18
- Step 4. Repair the ends 19
- Step 5. Purify the repaired DNA with the QIAquick PCR Purification Kit 20
- Step 6. Add 'A' Bases to the 3' end of the DNA fragments 21
- Step 7. Purify the sample with a Qiagen MinElute PCR Purification Column 22
- Step 8. Ligate the adapter 23
- Step 9. Purify the sample with the QIAquick PCR Purification Kit 24
- Step 10. Select size using a gel 25
- Step 11. Purify the gel 26
- Step 12. Amplify prepped library 28
- Step 13. Purify the sample 29
- Step 14. Check gDNA quantity and quality 30
- Step 15. Assess quality with Agilent 2100 Bioanalyzer 31

This section contains instructions for prepped library production specific to the Illumina single-read sequencing platform. It is intended for use with the Illumina prep kit (p/n FC-102-1001). The steps in this section differ from the Illumina protocol in the shear size, the use of the Covaris sample preparation system for gDNA shearing, and the gel purification. Other methods of gDNA shearing have not been validated.

Refer to the Illumina protocol *Preparing Samples for Sequencing Genomic DNA* (p/n 11251892 Rev. A) for more information.

NOTE

Make sure genomic DNA samples are of high quality with an OD 260/280 ratio ranging from 1.8 to 2.0.



2 Sample Preparation




Figure 1 Overall sequencing sample preparation workflow.

Table 7 Overview and time requirements

Step	Time
Illumina Prepped library Production	2 days
Library Hybridization	25 hours (optional 72 hours)
Bead preparation	30 minutes
Capture Selection and Washing	2 hours
DNA purification	30 minutes
Post-Hybridization Amplification	1 hour
PCR purification	30 minutes
Nanodrop and Bioanalyzer QC	1 hour

Step 1. Shear DNA

- 1 For every 3 µg DNA, add enough 1X Low TE Buffer in a LoBind tube for a total volume of 100 µL to dilute the genomic DNA.
- 2 Put a Covaris microTube into the loading station.
Keep the cap on the tube.
- 3 Set up the Covaris system.
Refer to the Covaris instrument user guide.
- 4 Use a tapered pipette tip to slowly transfer the 100 µL DNA sample through the pre-split septa.
Be careful not to introduce a bubble into the bottom of the tube.
- 5 Shear with the settings in [Table 8](#). The target peak for base pair size is 200.

Table 8 Covaris shear settings

Setting	Value
Duty Cycle	10%
Intensity	5
Cycles per Burst	200
Time	180 seconds
Set Mode	Frequency sweeping
Temperature	4°C

Step 2. Purify the sample with the QIAquick PCR Purification Kit

- 1 Add 500 μL of PB to the sample and mix well by pipetting.
- 2 Place a QIAquick spin column in a 2 mL collection tube.
- 3 Transfer the 600 μL sample to the QIAquick column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 4 Add 750 μL of buffer PE to the column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 5 Place the QIAquick column back in the 2 mL collection tube and spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 6 Transfer the QIAquick column to a new 1.5 mL collection tube to elute the cleaned sample. Add 30 μL of buffer EB (10 mM Tris·Cl, pH 8.5) directly onto the QIAquick filter membrane. Wait 60 seconds, then centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 7 Collect the eluate.

2 Sample Preparation

Step 3. Assess quality with Agilent 2100 Bioanalyzer

Step 3. Assess quality with Agilent 2100 Bioanalyzer

Use a Bioanalyzer DNA 1000 chip and reagent kit.

- 1 Check that the 2100 Bioanalyzer electrodes have been cleaned as instructed in the reagent kit guide.
- 2 Open the Agilent 2100 expert software (version B.02.02 or higher), turn on the 2100 Bioanalyzer and check communication.
- 3 Prepare the chip, samples and ladder as instructed in the reagent kit guide.
- 4 Load the prepared chip into the 2100 Bioanalyzer and start the run within five minutes after preparation.
- 5 Within the instrument context, choose the appropriate assay from the drop down list.
- 6 Start the run. Enter sample names and comments in the Data and Assay context.
- 7 Verify the results.




Figure 2 Analysis of sheared DNA using a DNA 1000 Bioanalyzer assay. The electropherogram shows a distribution with a peak size of 200 nucleotides.

Step 4. Repair the ends

- 1 In PCR tubes, strip tubes, or plates, prepare the reaction mix in [Table 9](#) for each prepped library, on ice. Mix well by gently pipetting up and down.

NOTE

If a precipitate forms in the T4 DNA ligase buffer with 10mM ATP, heat at 37°C for 5 minutes, and mix on a vortex mixer.

Table 9 End Repair

Reagent	Volume for 1 Library
DNA sample	29 µL
Nuclease-free water	46 µL
T4 DNA ligase buffer with 10mM ATP	10 µL
dNTP mix	4 µL
T4 DNA polymerase	5 µL
Klenow enzyme	1 µL
T4 PNK	5 µL
Total Volume	100 µL

- 2 Incubate in a thermal cycler for 30 minutes at 20°C.
 If you use a heated lid, make sure that the lid temperature does not exceed 50°C.

2 Sample Preparation

Step 5. Purify the repaired DNA with the QIAquick PCR Purification Kit

Step 5. Purify the repaired DNA with the QIAquick PCR Purification Kit

- 1 Add 500 μ L of PB to the sample and mix well by pipetting.
- 2 Place a QIAquick spin column in a 2 mL collection tube.
- 3 Transfer the 600 μ L sample to the QIAquick column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 4 Add 750 μ L of buffer PE to the column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 5 Place the QIAquick column back in the 2 mL collection tube and spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 6 Transfer the QIAquick column to a new 1.5 mL collection tube to elute the cleaned sample.
- 7 Add 32 μ L of buffer EB directly onto the QIAquick filter membrane. Wait 60 seconds, then spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 8 Collect the eluate.

Step 6. Add 'A' Bases to the 3' end of the DNA fragments

Step 6. Add 'A' Bases to the 3' end of the DNA fragments

- 1 In PCR tubes, strip tubes, or plates, prepare the reaction mix in [Table 10](#) for each prepped library, on ice. Mix well by gently pipetting up and down.

Table 10 Adding "A" Bases

Reagent	Volume for 1 Library
DNA sample	32 μ L
Klenow buffer	5 μ L
dATP	10 μ L
Klenow exo (3' to 5' exo minus)	3 μ L
Total Volume	50 μL

- 2 Incubate in a thermal cycler for 30 minutes at 37°C.

2 Sample Preparation

Step 7. Purify the sample with a Qiagen MinElute PCR Purification Column

Step 7. Purify the sample with a Qiagen MinElute PCR Purification Column

- 1 Allow the MinElute columns (stored at 4°C) to come to room temperature.
- 2 Add 250 µL of PB to the sample and mix well by pipetting.
- 3 Place a MinElute spin column in a 2 mL collection tube.
- 4 Transfer the 300 µL sample to the MinElute column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 5 Add 750 µL of buffer PE to the column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 6 Place the MinElute column back in the 2 mL collection tube and spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 7 Transfer the MinElute column to a new 1.5 mL collection tube to elute the cleaned sample. Add 10 µL buffer EB directly onto the MinElute filter membrane. Wait 60 seconds, then spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 8 Collect the eluate, which can be stored at 4°C.

Step 8. Ligate the adapter

This step uses a 10:1 molar ratio of adapter to genomic DNA insert, based on a starting quantity of 3 µg of DNA before fragmentation.

- 1 Prepare the reaction mix in [Table 11](#) on ice. Mix well by gently pipetting up and down.

Table 11 Ligation master mix

Reagent	Volume for 1 Library
DNA sample	10 µL
Nuclease-free water	4 µL
DNA ligase buffer	25 µL
Genomic Adapter oligo mix	6 µL
DNA ligase	5 µL
Total Volume	50 µL

- 2 Incubate for 15 minutes at room temperature.

2 Sample Preparation

Step 9. Purify the sample with the QIAquick PCR Purification Kit

Step 9. Purify the sample with the QIAquick PCR Purification Kit

- 1 Add 250 μL of PB to the sample and mix well by pipetting.
- 2 Place a QIAquick spin column in a 2 mL collection tube.
- 3 Transfer the 300 μL sample to the QIAquick column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 4 Add 750 μL of buffer PE to the column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 5 Place the QIAquick column back in the 2 mL collection tube and spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 6 Transfer the QIAquick column to a new 1.5 mL collection tube to elute the cleaned sample. Add 30 μL of buffer EB directly onto the QIAquick filter membrane. Wait 60 seconds, then centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 7 Collect the eluate.

Step 10. Select size using a gel

WARNING

Prolonged exposure to UV light can damage your DNA.

Excise as narrow a band as possible from the gel during gel purification.

CAUTION

Do not purify multiple samples on a single gel. Cross-contamination between libraries can occur.

- 1 Prepare a 4% NuSeive agarose gel with distilled water and TAE. Final concentration of TAE should be 1X. Use a large well comb to prevent overloading.
- 2 Add 3 μL of loading buffer to 8 μL of the low molecular weight DNA ladder.
- 3 Add 10 μL loading buffer to 30 μL of the DNA from the purified ligation reaction.
- 4 Load all of the ladder solution to one lane of the gel.
- 5 Load the entire sample in another lane of the gel, leaving at least a gap of one empty lane between ladder and sample. Electrophoresis parameters may need to be optimized.
- 6 Run the gel at 25 V for 17 hours.
- 7 Incubate the gel in SYBR Gold solution for 60 minutes.
- 8 View the gel on a Dark Reader transilluminator, which is a safer alternative to a UV transilluminator.
- 9 Use a clean scalpel or razor blade to excise the region of gel that contains DNA fragments in the 200 to 300 bp range.
Make sure you trim the gel slice to avoid the need to use two purification columns.

Step 11. Purify the gel

Use a Qiagen Gel Extraction Kit (Qiagen, p/n 28704) to purify the DNA from the agarose slices.

- 1** Weigh the gel slice. If the gel slice exceeds 400 mg, use 2 Qiagen QIAquick spin columns for purification.
- 2** Add 6 volumes of Buffer QG to 1 volume of gel (100 mg = 100 μ L).
You may need to do this in a 15 mL conical tube or in two tubes.
- 3** Incubate at room temperature for 10 minutes (or until the gel slice has completely dissolved). To help dissolve gel, mix the tube in a vortex mixer every 2 to 3 minutes during the incubation.
- 4** After the gel slice has dissolved completely, check that the color of the mixture is yellow (similar to Buffer QG without dissolved agarose).
- 5** Add 1 gel volume of isopropanol to the sample and mix.
For example, if the agarose gel slice is 100 mg, add 100 μ L isopropanol.
Do not spin the sample in a centrifuge at this time.
- 6** Place a QIAquick spin column in a provided 2 mL collection tube.
- 7** To bind DNA, apply the sample to the QIAquick column, and spin in a centrifuge for 1 minute at 17,900 x g (13,000 rpm).
The maximum volume of the column reservoir is 800 μ L. For sample volumes of more than 800 μ L, simply load and spin again.
- 8** Discard flow-through and place QIAquick column back in the same collection tube.
- 9** To wash, add 750 μ L of Buffer PE to QIAquick column and spin in a centrifuge for 2 minutes at 17,900 x g (13,000 rpm).
- 10** Discard the flow-through and spin the QIAquick column in a centrifuge for an additional 1 minute at 17,900 x g (13,000 rpm).
Residual ethanol from Buffer PE will not be completely removed unless the flow-through is discarded before this additional centrifugation.
- 11** Place QIAquick column into a clean 1.5 mL microcentrifuge tube.

12 To elute DNA, add 30 μL of Buffer EB to the center of the QIAquick membrane. Wait 60 seconds, then spin the column in a centrifuge for 1 minute.

If 2 Qiagen MinElute columns were used for a single sample, elute each one in 15 μL of EB, for a total of 30 μL per DNA sample.

Step 12. Amplify prepped library

This step uses PCR to selectively enrich those DNA fragments that have adapter molecules on both ends, and to amplify the amount of DNA in the library. The PCR is done with two primers that anneal to the ends of the adapters. Fourteen cycles of PCR are used.

- 1 Prepare the PCR reaction mix in [Table 12](#), on ice. Mix well by gently pipetting up and down.

Table 12 PCR Components

Reagent	Volume for 1 Library
DNA	1 μ L
Nuclease-free Water	22 μ L
PCR primer 1.1	1 μ L
PCR primer 2.1	1 μ L
DNA Polymerase (from Illumina Kit)	25 μ L
Total Volume	50 μL

- 2 Amplify using the following PCR program:

Table 13 PCR protocol

Step	Temperature	Time
Step 1	98°C	30 seconds
Step 2	98°C	10 seconds
Step 3	65°C	30 seconds
Step 4	72°C	30 seconds
Step 5		Repeat Step 2 through Step 4 for a total of 14 times.
Step 6	72°C	5 minutes
Step 7	4°C	Hold

Step 13. Purify the sample

- 1 Add 250 μL of PB to the sample and mix well by pipetting.
- 2 Place a QIAquick spin column in a 2 mL collection tub.
- 3 Transfer the 300 μL sample to the QIAquick column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 4 Add 750 μL of buffer PE to the column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 5 Place the QIAquick column back in the 2 mL collection tube and spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 6 Transfer the QIAquick column to a new 1.5 mL collection tube to elute the cleaned sample. Add 50 μL of buffer EB directly onto the QIAquick filter membrane. Wait 60 seconds, then spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 7 Collect the eluate, which can be stored at 4°C.

2 Sample Preparation

Step 14. Check gDNA quantity and quality

Step 14. Check gDNA quantity and quality

Use the NanoDrop ND-1000 UV-VIS Spectrophotometer (or equivalent) to assess gDNA concentration and purity.

- 1** Select **Nucleic Acid Measurement**, then select **Sample Type** to be **DNA- 50**.
- 2** Use 1.5 μL of EB buffer to blank the instrument.
- 3** Use 1.5 μL of each gDNA sample to measure DNA concentration. Record the gDNA concentration ($\text{ng}/\mu\text{L}$) for each sample. Calculate the yield (μg) by multiplying DNA concentration ($\text{ng}/\mu\text{L}$) by the sample volume and dividing by 1000.
- 4** Record the **A260/A280** and **A260/A230** ratios. High-quality gDNA samples should have an **A260/A280** ratio of 1.8 to 2.0, indicating the absence of contaminating proteins, and an **A260/A230** ratio of >2.0 , indicating the absence of other organic compounds such as guanidinium isothiocyanate, alcohol and phenol as well as cellular contaminants such as carbohydrates.

NOTE

A minimum of 500 ng of library is required for hybridization.

Step 15. Assess quality with Agilent 2100 Bioanalyzer

Use a Bioanalyzer DNA 1000 chip and reagent kit to assess the quality and size distribution of the PCR products.

- 1 Check that the 2100 Bioanalyzer electrodes have been cleaned as instructed in the reagent kit guide.
- 2 Open the Agilent 2100 expert software (version B.02.02 or higher), turn on the 2100 Bioanalyzer and check communication.
- 3 Prepare the chip, samples and ladder as instructed in the reagent kit guide.
- 4 Load the prepared chip into the 2100 Bioanalyzer and start the run within five minutes after preparation.
- 5 Within the instrument context, choose the appropriate assay from the drop down list.
- 6 Start the run. Enter sample names and comments in the Data and Assay context.
- 7 Verify the results.





Figure 3 Analysis of amplified prepped library DNA using a DNA 1000 Bioanalyzer assay. The electropherogram shows single peak in the size range of 200 to 300 nucleotides.

2 Sample Preparation

Step 15. Assess quality with Agilent 2100 Bioanalyzer



3 Hybridization

- Step 1. Hybridize the library 36
- Step 2. Prepare magnetic beads 40
- Step 3. Select hybrid capture with SureSelect 41
- Step 4. Desalt capture solution 42

In this step, you combine the prepped library with the blocking agents and the SureSelect Oligo Capture Library.

CAUTION

The ratio of SureSelect Oligo Capture Library to prepped library is critical for successful capture.






Figure 4 SureSelect Target Enrichment System Capture Process

Table 14 SureSelect Reagent Kit Components

Kit Component	250 RXN Kit	50 RXN Kit	Storage
SureSelect Hyb # 1	bottle	tube with orange cap	Room Temperature
SureSelect Hyb # 2	tube with red cap	tube with red cap	Room Temperature
SureSelect Hyb # 4	bottle	tube with black cap	Room Temperature
3M Sodium Acetate	tube with clear cap	tube with clear cap	Room Temperature
SureSelect Binding buffer	bottle	bottle	Room Temperature
SureSelect Wash Buffer #1	bottle	bottle	Room Temperature
SureSelect Wash Buffer #2	bottle	bottle	Room Temperature
SureSelect Elution Buffer	bottle	bottle	Room Temperature
SureSelect Neutralization Buffer	bottle	bottle	Room Temperature
SureSelect Hyb # 3	tube with yellow cap	tube with yellow cap	-20°C
SureSelect Block #1	tube with green cap	tube with green cap	-20°C
SureSelect Block #2	tube with blue cap	tube with blue cap	-20°C
SureSelect Block #3	tube with brown cap	tube with brown cap	-20°C
Rnase Block	tube with purple cap	tube with purple cap	-20°C
SureSelect GA PCR Primers	tube with clear cap	tube with clear cap	-20°C

CAUTION

You must avoid evaporation from the small volumes of the capture during the 24 hour or greater incubation.

If you want to use a different combination of thermocycler, lid temperature, plates or strips, and sealing method (strip caps or sealing tape), first test the conditions. Incubate 27 μ L of SureSelect Hybridization Buffers (without DNA) at 65°C for 24 hours (or longer, if applicable) as a test. Include buffer in each well that you might use, including those in the center and those on the edges. Check that you do not get extensive evaporation. Evaporation should not exceed 3 to 4 μ L.

For a partial list of tested options showing minimal evaporation, refer to [“Alternative Capture Equipment Combinations”](#) on page 50.

3 Hybridization

Step 1. Hybridize the library

Step 1. Hybridize the library

- 1 If prepped library concentration is below 147 ng/ μ L, use a vacuum concentrator to concentrate the sample.
- 2 Mix the components in Table 15 at room temperature to prepare the hybridization buffer.

Table 15 Hybridization Buffer

Reagent	Volume for 1 capture (μ L)	Volume for 5 captures (μ L)	Volume for 12 captures (μ L)
SureSelect Hyb # 1	25	125	275
SureSelect Hyb # 2 (red cap)	1	5	11
SureSelect Hyb # 3 (yellow cap)	10	50	110
SureSelect Hyb # 4	13	65	143
Total	49	245	539

- 3 If precipitate forms, warm the hybridization buffer at 65°C for 5 minutes.
- 4 Load 40 μ L of hybridization buffer per well into the “A” row of the PCR plate.

The number of wells filled in Row A is the number of libraries prepared. This example is for 12 captures.





Figure 5 Hybridization Buffer shown in Blue

Step 1. Hybridize the library

- 5** In a PCR plate, prepare prepped library for target enrichment:
- Adjust prepped library concentration to 500 ng in 3.4 μ L and add this to the “B” row in the PCR plate. Place each sample into a separate well.
 - Add 2.5 μ L of SureSelect Block #1 (green cap) to row B.
 - Add 2.5 μ L of SureSelect Block #2 (blue cap) to row B.
 - Add 0.6 μ L of SureSelect Block #3 (brown cap) to row B.
 - Mix by pipetting.
 - Seal the wells of rows “A” and “B” with caps and place the PCR plate in the thermocycler.
 - Run the following thermocycler program in [Table 16](#).
 - Incubate the prepped library and blockers at 95°C for 5 minutes.
 - Use a heated lid on the thermocycler at 105°C to hold the temperature of the plate on the thermocycler at 65°C.
- Make sure that the plate is at 65°C for a minimum of 5 minutes before you get to [step 7](#).

Table 16 PCR program

Step	Temperature	Time
Step 1	95°C	5 minutes
Step 2	65°C	∞

**Figure 6** Prepped library shown in red

3 Hybridization

Step 1. Hybridize the library

- 6** In a separate PCR plate, strip tubes, or tubes, prepare SureSelect Oligo Capture Library Mix for target enrichment:
 - a** Add 5 μ L of SureSelect Oligo Capture Library.
 - b** Add 1 μ L of nuclease-free water to the Capture Library.
 - c** Use nuclease-free water to prepare a 1:1 dilution of the RNase Block (purple cap).
 - d** Add 1 μ L of diluted RNase Block to each Capture Library, and mix by pipetting.
 - e** Add the Capture Library mix to the “C” row in the PCR plate.
 - f** For multiple samples, use a multi-channel pipette to load the Capture Library samples into the “C” row in the PCR plate (see [Figure 7](#) for positions).

Keep the plate at 65°C during this time.
 - g** Seal the wells with strip caps. Use a capping tool to make sure the fit is tight.
 - h** Incubate the samples at 65°C for 2 minutes.

Step 1. Hybridize the library

- 7** Maintain the plate at 65°C while you use a multi-channel pipette to take 13 μL of Hybridization Buffer from the “A” row and add it to the SureSelect Capture Library Mix contained in row “C” of the PCR plate for each sample.




Figure 7 SureSelect Capture Library, or “Baits”, shown in Green

- 8** Maintain the plate at 65°C while you use a multi-channel pipette to add 7 μL of each prepped library mix in row “B” to the hybridization solution in row “C”. Mix well by pipetting up and down 8 to 10 times.
- 9** Seal the wells with strip caps. Make sure all wells are completely sealed.
The hybridization mixture is now 27 μL .
- 10** Incubate the hybridization mixture for 24 hours at 65°C with a heated lid at 105°C.
Samples may be hybridized for up to 72 hours, but when you hybridize at longer periods, check that there is no extensive evaporation.

3 Hybridization

Step 2. Prepare magnetic beads

Step 2. Prepare magnetic beads

- 1 Prewarm SureSelect Wash Buffer #2 at 65°C in a circulating water bath for use in “Step 3. Select hybrid capture with SureSelect”.
- 2 Vigorously resuspend the Dynal (Invitrogen) magnetic beads on a vortex mixer. Dynal beads settle during storage.
- 3 For each hybridization, add 50 µL Dynal magnetic beads to a 1.5 mL microfuge tube.
- 4 Wash the beads:
 - a Add 200 µL SureSelect Binding buffer.
 - b Mix the beads on a vortex mixer for 5 seconds.
 - c Put tubes into a magnetic device, such as the Dynal magnetic separator (Invitrogen).
 - d Remove and discard the supernatant.
 - e Repeat [step a](#) through [step d](#) for a total of 3 washes.
- 5 Resuspend the beads in 200 µL of SureSelect Binding buffer.

Step 3. Select hybrid capture with SureSelect

- 1 Add the hybridization mixture directly from the thermocycler to the bead solution, and invert the tube to mix 3 to 5 times.
- 2 Incubate the hybrid-capture/bead solution on a Nutator for 30 minutes at room temperature.
Make sure the sample is properly mixing in the tube.
- 3 Separate the beads and buffer on a Dynal magnetic separator and remove the supernatant.
- 4 Resuspend the beads in 500 μ L SureSelect Wash Buffer #1 by mixing on a vortex mixer for 5 seconds.
- 5 Incubate the samples for 15 minutes at room temperature.
- 6 Wash the beads:
 - a Separate the beads and buffer on a Dynal magnetic separator and remove the supernatant.
 - b Mix the beads in prewarmed 500 μ L SureSelect Wash Buffer #2 on a vortex mixer for 5 seconds to resuspend the beads.
 - c Incubate the samples for 10 minutes at 65°C.
 - d Invert the tube to mix. (The beads may have settled.)
 - e Repeat entire [step a](#) through [step d](#) for a total of 3 washes.
Make sure all of the wash buffer has been removed.
- 7 Mix the beads in 50 μ L SureSelect Elution Buffer on a vortex mixer for 5 seconds to resuspend the beads.
- 8 Incubate the samples for 10 minutes at room temperature.
- 9 Separate the beads and buffer on a Dynal magnetic separator.
- 10 Use a pipette to move the supernatant to a new 1.5 mL microcentrifuge tube.
- 11 Add 50 μ L of SureSelect Neutralization Buffer.


Step 4. Desalt capture solution

In this step you desalt the capture solution with a Qiagen minElute PCR purification column.

- 1** Allow the MinElute columns (stored at 4°C) to come to room temperature.
- 2** Add 500 µL of PB to the sample and mix well by pipetting.
- 3** If pH indicator I was added to the PB buffer, check for the yellow color to make sure buffer PB pH is correct.

For more information on how to check buffer pH, refer to the Qiagen MinElute Handbook to check for proper pH. If needed, use 5 µL of the 3M Sodium Acetate (included in the kit) to adjust the pH to the proper range.

- 4** Place a MinElute spin column in a 2 mL collection tube.
- 5** Transfer the 600 µL sample to the MinElute column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 6** Add 750 µL of buffer PE to the column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 7** Place the MinElute column back in the 2 mL collection tube and spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 8** Transfer the MinElute column to a new 1.5 mL collection tube to elute the cleaned sample. Add 15 µL buffer EB directly onto the MinElute filter membrane. Wait 60 seconds, then spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 9** Collect the eluate (captured library), which can be stored at -20°C.



4 Post-Hybridization Amplification

Step 1. Amplify the sample 44

Step 2. Purify PCR 46

Step 3. Assess quality with Agilent 2100 Bioanalyzer 47

This chapter describes the steps to amplify, purify, and assess quality of the sample library.



4 Post-Hybridization Amplification

Step 1. Amplify the sample

Step 1. Amplify the sample

CAUTION

Do not use amplification enzymes other than Herculase II Fusion DNA Polymerase. Other enzymes have not been validated.

- 1 For each hybrid capture, prepare 2 amplification reactions. Mix the components in [Table 17](#) on ice, to create the amplification reactions. Mix well by gently pipetting up and down.

Table 17 Herculase Master Mix

Reagent	Volume for 2 reactions (including excess)	Volume for 10 reactions (including excess)	Volume for 24 reactions (including excess)
Nuclease-free water	91.25 μ L	401.5 μ L	912.5 μ L
5X Herculase II Reaction Buffer	25 μ L	110 μ L	250 μ L
dNTP mix (25 mM each)	1.25 μ L	5.5 μ L	12.5 μ L
SureSelect GA PCR Primers	2.5 μ L	11 μ L	25 μ L
Herculase II Fusion DNA Polymerase	2.5 μ L	11 μ L	25 μ L
Total	122.5 μL	539 μL	1,225 μL

- 2 Add 49 μ L of the mixture into a PCR plate or tubes for each capture.
- 3 Add 1 μ L of the Captured library to each Master Mix.

4 Place the tubes in a thermal cycler and run the program in [Table 18](#).

Table 18 PCR program

Step	Temperature	Time
Step 1	98°C	30 seconds
Step 2	98°C	10 seconds
Step 3	57°C	30 seconds
Step 4	72°C	30 seconds
Step 5		Repeat Step 2 through Step 4 for a total of 18 times.
Step 6	72°C	7 minutes
Step 7	4°C	Hold

4 Post-Hybridization Amplification

Step 2. Purify PCR

Step 2. Purify PCR

- 1 Combine duplicate amplifications for a total of 100 μL .
- 2 Add 500 μL of PB to the sample and mix well by pipetting.
- 3 Place a QIAquick spin column in a 2 mL collection tube.
- 4 Transfer the 600 μL sample to the QIAquick column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 5 Add 750 μL of buffer PE to the column. Spin the sample in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm). Discard the flow-through.
- 6 Place the QIAquick column back in the 2 mL collection tube and spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 7 Transfer the QIAquick column to a new 1.5 mL collection tube to elute the cleaned sample. Add 30 μL of buffer EB directly onto the QIAquick filter membrane. Wait 60 seconds, then spin in a centrifuge for 60 seconds at 17,900 x g (13,000 rpm).
- 8 Collect the eluate, which can be stored at -20°C .
- 9 Measure DNA concentration with a spectrophotometer (NanoDrop Technologies).

Step 3. Assess quality with Agilent 2100 Bioanalyzer

Use a Bioanalyzer DNA 1000 chip and reagent kit.

- 1 Check that the 2100 Bioanalyzer electrodes have been cleaned as instructed in the reagent kit guide.
- 2 Open the Agilent 2100 expert software (version B.02.02 or higher), turn on the 2100 Bioanalyzer and check communication.
- 3 Prepare the chip, samples and ladder as instructed in the reagent kit guide.
- 4 Load the prepared chip into the 2100 Bioanalyzer and start the run within five minutes after preparation.
- 5 Within the instrument context, choose the appropriate assay from the drop down list.
- 6 Start the run. Enter sample names and comments in the Data and Assay context.
- 7 Verify the results.
- 8 Continue to sequencing.




Figure 8 Analysis of Amplified Capture DNA using a DNA 1000 Bioanalyzer assay. The electropherogram shows a single peak in the size range of 200 to 300 nucleotides.

4 Post-Hybridization Amplification

Step 3. Assess quality with Agilent 2100 Bioanalyzer



SureSelect Target Enrichment
Protocol

5 Reference

Alternative Capture Equipment Combinations 50

This chapter contains reference information.



Alternative Capture Equipment Combinations

Table 19 lists combinations of thermocycler, lid temperature, plates or strips, and sealing method (strip caps or sealing tape) other than those used in this protocol that have shown minimal evaporation.

Refer to this list for additional of equipment combination options for hybridization. Note that minimal evaporation is needed to ensure good capture results.

Table 19 Tested options that show minimal evaporation

PCR Machine	Plate/ Strips	Cover	Comments
Stratagene Mx3005P QPCR	Mx3000P Strip Tubes (401428)	MX3000P Optical Strip Caps (401425)	Heated Lid
Stratagene Mx3005P QPCR	MicroAmp Optical 96-well reaction plate (N801-0560)	MicroAmp Clear Adhesive Film (4306311)	Heated Lid; ABI compression pad (4312639) Use two layers of film.
ABI GeneAmp 9700	MicroAmp Optical 96-well Reaction Plate (N801-0560)	MicroAmp Caps (8caps/strip) (N801-0535)	Heated Lid
ABI Veriti (4375786)	MicroAmp Optical 96-well Reaction Plate (N801-0560)	MicroAmp Clear Adhesive Film (4306311)	Heated Lid; ABI compression pad (4312639) Use two layers of film.
Eppendorf Mastercycler	Eppendorf 8-Tube PCR Tubes	Attached lids	Lid heating set to 75C
ABI Veriti (4375786)	Stratagene strip tubes 410022 (Mx4000)	Stratagene Strip cap domed 410096 (Robocycler)	Heated Lid
ABI Veriti (4375786)	Stratagene strip tubes 410022 (Mx4000)	Stratagene Optical cap 401425 (Mx3000/3005)	Heated Lid
MJ Research PTC-200	Stratagene strip tubes 410022 (Mx4000)	Stratagene Optical cap 410024 (Mx4000)	Heated Lid
MJ Research PTC-200	Stratagene strip tubes 410022 (Mx4000)	Stratagene Optical cap 401425 (Mx3000/3005)	Heated Lid

Table 19 Tested options that show minimal evaporation (continued)

PCR Machine	Plate/ Strips	Cover	Comments
MJ Research PTC-200	Stratagene 96-well Plate 410088 (Mx3000/3005)	Stratagene Optical cap 401425 (Mx3000/3005)	Heated Lid
MJ Research PTC-200	Stratagene 96-well Plate 410088 (Mx3000/3005)	Stratagene Plate sealers 400774-15	Heated Lid 2 layers of plate sealer

5 Reference
Alternative Capture Equipment Combinations

www.agilent.com

In This Book

This guide contains information to run the SureSelect Target Enrichment protocol.

© Agilent Technologies, Inc. 2009

Printed in USA
Version 1.2, April 2009



G3360-90010



Agilent Technologies