

PROTOCOL

Preparation of Felodipine Nanosuspension (2 ml) by crystallization through AFA technology

Summary of Operating Conditions

Table 1. Summary of Operating Conditions

Peak Incident Power	175 W
Duty Factor	20%
Cycles per Burst	1000
Duration	1200 seconds
Water Bath Temperature	18°C
Power Mode	Frequency sweeping
Degassing Mode	Continuous
Volume	2 mL

Supplies

Table 2. Equipment List

Part Name	Description	Part Number
Sample Vessel	Tube & Cap 12x24mm	520056
Sample Holder	Holder 12x24 Tube	500199
Focused Ultrasonicator	Covaris S220x	S220x
Sample	Felodipine	Sigma, F9677
	Dimethylacetamide	Sigma, 38840
Vehicle	0.2% Polyvinylpyrrolidone (PVP30)	BASF
	0.25 mM Sodium Lauryl Sulfate (SLS)	Sigma

Operating Conditions

1. Fill the tank with fresh deionized water to proper fill line. The S220x should be equipped with a graduated water level label. If the tank lacks this label, please contact Covaris. During treatments, the tube should be partially immersed in the water to ensure a good acoustic path from the AFA transducer.
2. Degas water for the recommended 30 minutes or more. To maintain degassed water, keep the pump continuously on during operation and sample processing. Do not turn the pump off.
3. Set the chiller to the proper temperature, as listed for “Water Bath Temperature” in Table 1.

Recommendations Specific for Nanosuspension Formation

The Covaris AFA process is highly reproducible, however steps should be taken to ensure the best results. The bath water is employed to couple acoustic energy to the sample vessel, thus attention must be paid to the following water treatment attributes to obtain the best results:

1. *Purity*: When applying acoustics in rate-limited applications, foreign materials such as algae and particulates may scatter the high frequency focused acoustic beam. Bath water should be pure distilled or DI water, changed daily or cleansed by a Covaris Water Conditioning System.
2. *Degas Level*: Similarly, insufficient degas levels within the bath may result in poor acoustic coupling. System degas pumps should be run in advance of and during AFA treatments, as detailed in instrument User Manuals.
3. *Temperature*: Warmer temperatures promote less forceful collapse of acoustic cavities within the sample fluid. Bath temperature (as reported by SonoLAB software) should therefore be closely controlled and matched run-to-run and day-to-day. Employ the temperature alert feature in SonoLAB to warn of a failure to maintain control of bath temperature.
4. *Level*: Attention should be paid to maintaining a consistent water level, according to published protocols. If using a Covaris Water Conditioning System, check levels daily to restore water lost to evaporation.

In summary, when employing the Covaris AFA, control and verification of treatment attributes and water quality will reduce variance and promote consistent, satisfactory results.

Method

1. Set up the Covaris S220x at the appropriate temperature following the operating conditions above.
2. Add 20 μL of a 100 mM stock solution of Felodipine in dimethylacetamide to the sample vessel.
3. Add 1980 μL of the vehicle (PVP30 and SLS) to the sample vessel, and cap the tube. This will result in a 99:1 ratio of the vehicle to the drug solution (i.e. the drug concentration will be 1 mM). Be careful not to introduce a bubble into the bottom of the tube. This may happen if the water is added too quickly.
CAUTION: The bottom of the tube is in the acoustic field. Therefore, a bubble in the sample will deflect energy and induce variable results.
4. Carefully load the sample vessel into the appropriate holder, and insert the holder into the S220x instrument.
5. Initiate and Run process according to the operating conditions specified in Table 1.

PROTOCOL

Description

A Malvern Zetasizer – 90 NS-90 may be used to analyze the nanosuspension formed using the Covaris S220x. The sample prepared above will be utilized for analysis.

Supplies

Table 3. Supplementary Equipment List


Part Name	Part Number
Malvern Zetasizer ZS-90	ZEN3690
Sample Cuvettes	DTS0012

Method

1. Set up the Malvern Zetasizer according to its setup instructions.
2. Add 1 mL of processed sample to the cuvette.
3. Cap the cuvette.
4. Shake the cuvette by hand until the sample is dispersed evenly in the cuvette without air bubbles.
5. Place in Zetasizer instrument and run Volume Distribution Analysis.

Typical Output Readings: Intensity

Figure 1. Size Distribution by Intensity



Z-average = 153.7 nm; polydispersity index = 0.107

PROTOCOL

Typical Output Readings: Variation

Table 4. Variation in Results

Process	PDI	Z-average (nm)
1	0.107	153.7
2	0.099	135.7
3	0.136	154.6
4	0.178	172.6
Average	0.130	154.2
SD	0.031	13.1
CV	24%	8%

Figure 2. Variation in Nanosuspension Formation Results

