

truXTRAC® FFPE Total Nucleic Acid Kits

Simultaneously Extract NGS-quality DNA and RNA from FFPE Tissue Samples


The truXTRAC Total Nucleic Acid kits use the patented Adaptive Focused Acoustics® (AFA®) technology to actively and concurrently extract both DNA and RNA from formalin-fixed, paraffin-embedded (FFPE) tissue in an aqueous buffer. In addition to automated magnetic bead-based formats, manual column-based purification formats are also available. These kits have been designed to maximize recovery of intact nucleic acids, allowing researchers to detect low frequency mutations and rare gene fusion events.

Benefits


- Higher DV₂₀₀ scores
- Lower QNS rates
- Superior coverage depth
- Improve gene fusion detection
- Deparaffinization and rehydration
- Automation-friendly
- DNA/RNA co-extraction
- Minimize sample input required
- Safe: No organic solvent
- Higher extraction yield

Applications

- Whole Genome Sequencing
- Whole Exome Sequencing
- Amplicon Sequencing
- RNA-seq
- qPCR
- Droplet digital PCR (ddPCR)



Simple, Standardized Workflow



By using the Covaris truXTRAC FFPE DNA and RNA kits, OmniSeq testing requires less tissue and results in fewer specimen failures when performing comprehensive genomics profiling on solid tumors.


RNA DV₂₀₀ Score Comparison

Figure 1. Illumina uses 30% DV₂₀₀ cut-off. DV₂₀₀: Distribution Value 200, which is the percentage of RNA fragments larger than 200 nucleotides. (Levy, S. The ORIEN Project- Large Scale Integration of FFPE Material into a Multisite Research Alliance. Oral presentation at AACR 2019; April 1, 2019)

Whole Genome Sequencing



Figure 2. Whole Genome Sequencing Results (Chromosome 19 Coverage) Comparison Covaris truXTRAC, Competitor and fresh frozen tissue. Sequencing track from QIAGEN (red), Covaris (green) extracted FFPE kidney DNA, and fresh frozen (indigo) samples were loaded on the IGV viewer, and the coverage analyzed for chromosome 19. Coverage of >10x are indicated in dark colors, coverage of <10X are indicated in light colors. Chromosomal view of coverage indicates that Covaris extracted DNA quality resembles that of the DNA extracted from fresh frozen tissue.

Ordering Information

Part Number	Product Name	Description
520246	truXTRAC FFPE total NA Kit - Magnetic Bead	Kits contains reagents and consumables required for the sequential extraction and purification of DNA and RNA from FFPE tissue samples. Optimized to work with standard FFPE inputs.
520255	truXTRAC FFPE total NA Kit Plus Kit - Magnetic Bead	
520220	truXTRAC FFPE total NA Kit - Column	Column-based reagent kits that contains the components needed to extract and purify DNA and RNA FFPE tissue samples. Optimized to work with standard FFPE inputs, this kit is provided with the microTUBE-130 for efficient paraffin emulsification and tissue rehydration.
520252	truXTRAC FFPE total NA Plus Kit - Column	