

DNA Extraction of Lung Cancer Samples for Advanced Diagnostic Testing

Conroy JM^{1,2}, Hanchun T. DeFedericis¹, Angela Stout¹, Kiersten Marie Miles¹, Burgher B¹, Papanicolau-Sengos A¹, Morrison CD^{1,2*}
OmniSeq, Inc¹, Buffalo, NY; Roswell Park Cancer Institute², Buffalo, NY *Carl.Morrison@OmniSeq.com

Introduction

It is apparent that genetic testing leads to improved survival rates for lung cancer patients. Patients whose tumor DNA harbor oncogenic driver mutations respond favorably to drug treatments specifically tailored to those drivers. Currently, these lung targeted therapies are limited to mutations in 15-20 genes, but the number will grow as evidence emerges supporting novel drug-mutation interactions. A major challenge today is detecting therapeutically actionable driver mutations in lung cancer biopsies, most of which are sub-optimally preserved and contain limited amounts of evaluable tumor tissue.

Background

Lung cancer biopsies are routinely formalin fixed and paraffin embedded (FFPE) for histological evaluation by anatomical pathologists. This process preserves the tumor morphology and cellular features required for proper staining and microscopic review. However, this practice presents numerous challenges for the extraction of high quality DNA for genomic testing. An extraction process that consistently produces sufficient DNA yield and fragment size from these difficult but most precious tissue samples is a requirement for any Molecular Pathology laboratory. The data presented here compares the quantity and quality of DNA extracted using two methods, QIAGEN and Covaris, and the success of downstream molecular diagnostic testing platforms routinely utilized for detecting somatic alterations including single nucleotide variants (SNV), insertion/deletions (indel), copy number (CN) amplification, microsatellite instability (MSI) and tumor mutational burden (Table 1).

Alteration	Gene	Method	Platform	DNA requirement (ng)
Mutation	BRAF	Sanger sequencing	ABI 3500	25-50
	EGFR	Pyrosequencing	PyroMark Q24	
	KRAS	single analyte		
Mutation/amplification/fusion	AKT1	Multi-analyte: NGS and digital detection	Illumina MiSeq Ion Torrent PGM, S5XL Nanostring nCounter	10-500
	ALK			
	BRAF			
	DDR2			
	EGFR			
	ERBB2			
	FGFR1			
	KRAS			
	MAP2K1			
	MET			
	NRAS			
	PIK3CA			
	PTEN			
RET				
ROS1				
Tumor mutational burden	> 400	Multi-analyte: NGS	Illumina MiSeq Ion Torrent S5XL	10-200

Table 1: Commonly tested alterations in lung cancer with method and DNA requirements.

Methods

FFPE tumor samples from a variety of tumor types, including lung, were macro-dissected using 14-gauge needles, with one core extracted using the Covaris truXTRAC FFPE DNA isolation method and the other matched core using the QIAGEN DNeasy tissue kit (Figure 1). All samples were processed using manufacturer's recommended instructions. DNA metrics were measured using Qubit and NanoDrop for yield and purity, followed by fragment size estimation on a 2100 BioAnalyzer.



A subset of matched DNA sample pairs were used as template for mutation detection using pyrosequencing (PyroMark), Ion AmpliSeq (PGM) and Illumina TruSeq Custom Amplicon (MiSeq) NGS, and for copy number using the NanoString nCounter system. An additional set of Covaris prepared DNA samples were tested at 1:2 and 1:4 dilutions. OmniSeq TargetSM for Lung, a NYS CLEP approved test, was used for all NGS and nCounter analyses.

Results

DNA yields and fragment lengths were substantially higher for Covaris extracted samples as compared to QIAGEN when measured by Qubit (Table 2) and Bioanalyzer electrophoresis (Figure 2). A higher degree of successful advanced molecular diagnostic test results was also observed for the truXTRAC DNA samples, especially for the MiSeq NGS system (improved clustering and coverage) and nCounter platform (improved counts) that prefer longer fragment lengths than PGM NGS.

Sample #	Disease Type	NanoDrop				Qubit	
		QIAGEN yield (ug)	Covaris yield (ug)	QIAGEN: DNA 260/280	Covaris: DNA 260/280	QIAGEN yield (ug)	Covaris yield (ug)
1	CRC	9.72	23.2	2.04	1.81	0.93	6.56
2	Melanoma	25.03	1.11	1.39	1.59	0.43	0.08
3	Breast CA	6.38	12.92	2.00	1.85	0.77	4.46
4	Breast CA	12.17	17.07	2.05	1.83	1.44	5.30
5	Breast CA	13.66	17.28	2.01	1.89	3.28	5.02
6	Melanoma	13.14	20.36	2.02	1.81	0.81	4.24
7	NSCLC	8.9	16.89	2.00	1.85	0.86	4.54
8	CRC	19.02	30.07	2.05	1.86	3.38	6.76
Average		13.50	17.36	1.95	1.81	1.49	4.62

Table 2: DNA yields of matched FFPE tissue core extractions (1mm).

Results




Figure 2: Tissue mass (# cores), DNA yield (ug), fragment traces, and NGS and nCounter CN test result success for six NSCLC samples extracted with QIAGEN (left) and Covaris (right) methodologies.

Covaris extracted DNA samples were also successfully validated for use with single analyte PyroMark assays, and when diluted 25-50% for NGS analysis (Table 3).

Sample #	PyroMark Test	QIAGEN result	Covaris result	Sample #	DNA Dilution	NGS QC	NGS SNV calls (concordance)
22	BRAF	mutation ND	mutation ND	27	1:1	pass	32
	EGFR	mutation ND	mutation ND		1:2	pass	32 (100%)
	KRAS	mutation ND	mutation ND		1:4	pass	32 (100%)
24	JAK2	mutation ND	mutation ND	28	1:1	pass	27
	BRAF	p.Val600Glu	p.Val600Glu		1:2	pass	27 (100%)
	EGFR	mutation ND	mutation ND		1:4	pass	27 (100%)
26	KRAS	mutation ND	mutation ND	29	1:1	pass	24
	JAK2	mutation ND	mutation ND		1:2	pass	24 (100%)
	BRAF	mutation ND	mutation ND		1:4	pass	24 (100%)
31	EGFR	mutation ND	mutation ND	31	1:1	pass	24 [#]
	KRAS	p.Gly12Asp	p.Gly12Asp		1:2	pass	25 (96%)
	JAK2	mutation ND	mutation ND		1:4	pass	25 (96%)

[#] RET c.2712C>G silent mut, VAF 0.40 not detected

Table 3: PyroMark results for QIAGEN and Covaris DNA (left) and OmniSeq Target NGS results for serially diluted Covaris DNA (right).

Conclusions

FFPE tumor samples prepared using the Covaris truXTRAC isolation kit provides an efficient system for generating high quality DNA samples even from lung cancer specimens that previously failed testing with QIAGEN extraction. The combination of improved yield and fragment size measured for nearly every sample tested suggests that even smaller biopsies can now be collected and extracted for advanced diagnostic testing. Additionally, DNA requirements for NGS based testing can be reduced 25-50% without loss of assay performance or analytic sensitivity as validated in our CLIA laboratory.