

15 μ l DNA Shearing with microTUBE-15™

NEW

microTUBE-15

- 15 μ l DNA shearing for NGS library prep
- Fast, simple workflow
- High DNA recovery
- Isothermal process eliminates heat-induced bias
- Accurate, tunable DNA fragmentation
- Robust and reproducible results day to day, sample to sample


8 microTUBE-15 Strip



For more details scan or visit:

[http://covarisinc.com/
products/afa-tubes-and-vials/
microtube-15/](http://covarisinc.com/products/afa-tubes-and-vials/microtube-15/)

Accurate: Tunable shearing length



ORDERING INFORMATION

Part #	Product Name	Description
520145	microTUBE-15 AFA Beads Screw-Cap (25)	<ul style="list-style-type: none"> • microTUBE-15 with AFA Beads™ for fully controllable and repeatable DNA shearing in very low sample volume. • 2D barcodes on each microTUBE-15 provide individual sample tracking.
520159	8 microTUBE-15 AFA Beads Strip V2 (12)	<ul style="list-style-type: none"> • Strip of microTUBE-15 with AFA Beads™ for fully controllable and repeatable DNA shearing in very low sample volume. • 2D barcodes on each strip provide sample tracking.

ACCESSORIES

520059	Centrifuge microTUBE Adapters (25)	Required to hold the microTUBE-15 in a standard microcentrifuge rotor designed for 1.5 or 2 mL microcentrifuge tubes. All microTUBE-15 protocols specify an initial centrifugation step.
500414	M220 Holder XTU	For use with a Covaris M220 Focused-ultrasonicator. Requires an Insert.
500420	M220 Holder XTU Insert microTUBE 15 µL	For use with a Covaris M220 Focused-ultrasonicator and Holder XTU.
500427	S-Series Holder Screw-Cap LV	For use with a Covaris S220 Focused-ultrasonicator.
500308	E-Series Rack 24 Place microTUBE Screw-Cap	For use with a Covaris E220 Focused-ultrasonicator.
500444	E-Series Rack 12 Place 8 microTUBE Strip V2	For use with a Covaris E220 Focused-ultrasonicator.

STARTER KITS

500417	M220 NGS microTUBE-15 Starter Kit	Includes: <ul style="list-style-type: none"> • M220 Holder XTU • M 220 Holder XTU Insert microTUBE 15 µL • microTUBE-15 AFA Beads Screw-Cap (25) • Centrifuge microTUBE Adapters (25) • microTUBE Prep Station Snap & Screw Cap
500418	S220 NGS microTUBE-15 Starter Kit	Includes: <ul style="list-style-type: none"> • S-Series Holder Screw-Cap LV • microTUBE-15 AFA Beads Screw-Cap (25) • Centrifuge microTUBE Adapters (25) • microTUBE Prep Station Snap & Screw Cap
500419	E220 NGS 24 microTUBE-15 Starter Kit	Includes: <ul style="list-style-type: none"> • Rack 24 Place microTUBE Screw-Cap • microTUBE-15 AFA Beads Screw-Cap (25) • Centrifuge microTUBE Adapters (25) • microTUBE Prep Station Snap & Screw Cap
500446	E220 NGS 96 microTUBE-15 Starter Kit	Includes: <ul style="list-style-type: none"> • E-Series Rack 12 Place 8 microTUBE Strip V2 • 8 microTUBE-15 AFA Beads Strip V2 (12)

Validated for M220, S220, and E220. For earlier models please contact Covaris.

USA: Covaris, Inc. • Tel: +1 781-932-3959 • Fax: +1 781-932-8705 • Email: customerservice@covarisinc.com • Web: www.covarisinc.com

EUROPE: Covaris Ltd. • Tel: +44 (0)845 872 0100 • Fax: +44 (0)845 384 9160 • Email: eucustomerservice@covarisinc.com • Web: www.covarisinc.com

PART NUMBER: M020027 REV A | EDITION APRIL 2016

INFORMATION SUBJECT TO CHANGE WITHOUT NOTICE | FOR RESEARCH USE ONLY | NOT FOR USE IN DIAGNOSTIC PROCEDURES | COPYRIGHT 2014 COVARIS, INC.